

BOTTLESHOWER GUIDE TO RUN TIMES

TYPE OF BOTTLE	LITRES	WITHOUT	WITH SHOWER HEAD**	WITH TAP HEAD**	NOTES
		Poured	Run time*	Run time*	
Life	500 ML	9 Sec	40 Sec	1 Min	Wide neck/ Flimsy Bottle
Wenlock Spring	500 ML	9 Sec	40 Sec	1 Min	Narrow Neck/Sturdy Bottle
Princess Gate	750 ML	11 Sec	48 Sec	1 Min 28 Sec	Narrow Neck/Sturdy Bottle
Evian	1 Litre	12 Sec.	1 Min. 15 Sec.	2 Min	Wide neck/ Sturdy Bottle
Buxton	1.5 Litre	16 Sec	1 Min 40 Sec	2 Min 33 Sec	Wide neck/ Sturdy Bottle
Highland Spring	2 Litre	-	2 Min	3 Min 20 Sec	Narrow neck/ Flimsy Bottle (Green)
Saskia	2 Litre	-	2 Min 10 sec	3 Min 15 Sec	Wide neck/ Sturdy Bottle
Celtic Spring	5 Litre	-	4 Min 30 Sec	8 Min 45 Sec	Large neck 39 - 42mm range
Re-usable Bottle	2.2 Litre	-	2 Min 15 Sec	3m 40 Sec	Large neck 39 - 42mm range ridged with handle
Specialist Bottles					
MOD/ Nato Osprey Pattern 58 bottle	1 Litre +	-	1 Min 20 Sec	2 Min 15 Sec	Large neck 39 - 42mm range
MOD/ Nato Water Carrier	20 Litre Capacity !0 Litre Test	-	10 Min	Not Relevant	Large neck 39 - 42mm range half capacity test.
Baxter Sterile water	1 Litre	13 Sec	1 Min	1 Min 50 Sec	Specialist bottle neck / Click fit
Fresenius Kabi Sterile water	1 Litre	13 Sec	1 Min 5 Sec	1 Min 55 Sec	Specialist bottle neck / Push fit
ALPHA XXL bottle	2 Litre	-	2 Min 15 sec	3 Min 10 Sec	Use with silicon bung, excellent robust bottle design for reuse
ALPHA XXL bottle	1.2 Litre	-	1 Min 20 Sec	2 Min 25 Sec	Use with silicon bung, excellent robust bottle design for reuse
			* Shower not dripping / average times		
			** Times will vary slightly under different conditions due to bottle shapes, these tests are not scientific but a fair report.		